

1. Part Number: A904-165

2. Photograph:



3. Product Information:

Description	Funnel with Strainer
Material of Construction:	304 stainless steel
Side Handle:	Spot welded to cone section of funnel
Surface Finish:	Better than 0.5 microns Ra
Top Diameter:	165mm
Outer Diameter of Spout:	12mm
Overall Height:	165mm
Diameter of Holes in strainer:	4mm
Is strainer removeable?	No
Nominal Weight:	212g
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 13-Oct-20



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.