


**1. Part Number:** 1220B-900

**2. Photo:**



**3. Product Information:**

Description	Pocket Sampler 900mm
Material of Construction:	316L stainless steel
External Polish:	Better than 1.0 microns Ra
Internal Pocket Finish:	Fine machined finish
Outer Diameter of Outer Tube:	25 mm
Overall Length:	900 mm
Pocket Length:	135 mm
Pocket Width:	11 mm
Nominal Pocket Volume:	20 ml
Number of Pockets	3
Pocket Positions:	



Nominal Weight of Sampler:	2.7 Kg
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 8-Feb-21



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.