

**1. Part Number:** 1430S-1500

**2. Photo:**



### 3. Product Information:

Description

Mono Thief 1500mm

Materials:

Outer Rod: 316L stainless steel  
Tip: 316L stainless steel  
Handle: 316L stainless steel  
Inner Rod: 316L stainless steel

Outer Diameter of Outer Tube:

25 mm

Overall Length:

1555 mm

Pocket Length:

135 mm

Pocket Width:

11 mm

Pocket Volume:

20 ml

Number of Pockets

1

Pocket Position:



External Finish:

<1 Micron Ra

Internal Pocket Finish:

Fine machined finish

Nominal Weight of Sampler:

5.3 Kg

BSE/TSE statement:

All polishing compounds used in the manufacture of this sampler are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Jan-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.