




1. Part Number: 2015A-10

2. Photo:



3. Product Information:

Description	Cohesive Sampler 10ml
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.8 microns Ra
Nominal Sample Volume:	10ml
Overall Length:	1030mm
Diameter of Outer Tube:	25mm
Scraper Length	110mm
Nominal Weight of Sampler:	1.9Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.