


1. Part Number: 2220B-600

2. Photo:



3. Product Information:

Description	Cohesive Pocket Sampler 600mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Overall Length:	600mm
Outer Diameter of Outer Tube:	25mm (plus 10mm approx. for scrapers)
Inner Diameter of Inner Tube:	17mm
Pocket Length:	115mm
Pocket Width:	11mm
Pocket Volume:	15ml
Number of Pockets	3
Pocket Positions:	



External Finish:	< 0.8 Microns Ra
Internal Pocket Finish:	Fine machined finish
Nominal Weight of Sampler:	1.7 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.