


1. Part Number: 5035T-1501

2. Photo:



3. Product Information:

Description

Liquid Master PTFE – 1500mm

Materials of Construction:

PTFE & FEP (for product contact parts)

- FDA 21 CFR 177.1550
- USP Class VI
- EU regulation 10/2011
- EC regulation 1935/2004

Outer Diameter of Outer Tube:

25mm

Internal Diameter of Outer Tube:

19mm

Max. Diameter of Tip:

32mm

Diameter of hole in end of tip:

10mm

Overall Length:

1500mm

BSE/TSE Free:

Yes

Nominal Weight:

860g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.