

 SAMPLING SYSTEMS	Product Specification	Specification Number	8095B-250S
		Revision	01

Part No. 8095B-250S

Description Disposable Powder Spatula Blue – 250mm

Material of Construction: Polypropylene (PP) – rigid, blue

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Method of Construction: Single piece, injection moulded

Nominal Length: 330mm

Blade Length: 246mm

Blade Width: 18mm

Nominal Weight: 22g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 21-Apr-17
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.