

AF 103 XSafe

Self Contained Ice Machine 120 kg

Scotsman®
Ice Systems

Item #:

Project:

Quantity:

		FEATURES <ul style="list-style-type: none">• Produces granular shape ice flakes, 75% hardness.• Electronically controlled operation, with self-diagnostic functions with external alarm lights providing operational information for early alert and fast diagnosis of operating issues.• External Components constructed from corrosion-resistant stainless steel crafted for optimal aesthetic appeal through superior fit and finish.• Ergonomically convenient, slanted & disappearing storage bin door.• Front access for easier maintenance.• Front-in, front-out airflow allows for built-in installations.• XSafe: natural sanitation system, integrated inside the ice maker and operating automatically 24/7 to keep equipment clean & safe.	
 <div>AF 103 Flake ice Residual water content 25%</div>	CONDENSING SYSTEM Air cooled Water cooled	INTERNAL BIN CAPACITY 30 kg	
	REFRIGERANT GAS R134a		
	VOLTAGE V/Hz/ph 230/50/1 230/60/1		
CERTIFICATIONS     		OPERATING REQUIREMENTS	
		Minimum	Maximum
Air temperature		10	43
Water temperature		5	38
Water pressure		1 bar (14 psi)	5 bar (70 psi)
Electrical voltage		-10%	+10%
IMPORTANT NOTICE: Models and specifications are subject to change without notice. This spec sheet is meant for commercial purpose only. For technical documentation please refer to our service manuals.			


Self Contained Ice Machine 120 kg

10 November 2020

AF 103 XSafe

Self Contained Ice Machine 120 kg

Item #:
Project:
Quantity:



UNIT DATA					SHIPPING DATA						
Size (W x D x H)		592 x 622 x 1078 mm			Carton (W x D x H)		670 x 690 x 1180 mm				
Net weight		66 kg			Weight		74 kg				
		24 h ice production kg °C Amb. / °C Water			Energy consumption (*)		Water usage (*)	Instant power	Compressor	Circuit wires	Max. fuse size
Version	Voltage	10°C/10°C	21°C/10°C	32°C/21°C	kWh/ 100 kg	kWh/ 24h	l/h	W	W	No. x Ømm²	A
		kg	kg	kg							
AF 103 AS	230/50/1	120	120	93	12.1	11.5	4.2	550	1200	3 x 1.5	10
AF 103 WS	230/50/1	110	108	90	11.8	10.5	32.0	440	1200	3 x 1.5	10
AF 103 AS	230/60/1	131	123	97	12.0	12.0	4.0	540	1400	3 x 1.5	10

(*) Data refer to 32°C Amb. / 21°C Water temperature conditions



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.