


Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Bibbypet Bottle top Dispensers BD10, BD25

- Economical dispensing of non-aggressive liquids
- Low-cost replaceable piston/cylinder and valves
- Easy cleaning and maintenance
- Simple to operate
- Automatic refill



The BibbyPet bottle top dispenser provides a low cost option for liquid dispensing without compromising on accuracy or reliability. The Bibbypet's innovative design utilises materials suitable for routinely dispensing many common laboratory liquids; dilute acids, buffers, biological solutions, polar solvents, even strong alkaline solutions.


* The values of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20°C) and with smooth operation.

Model	Volume, ml	Subdivision, ml	Accuracy, %	Coefficient of variation, %
BD10	1 to 10	0.2	1.2	0.2
BD25	2.5 to 25	0.5	1.2	0.2

Accessories

Polypropylene adapters

Model	Adapter thread, mm	For bottle thread, mm	Pack Qty
BD3222	GL32	GL22	1
BD3228	GL32	GL28	1
BD3230	GL32	GL30	1
BD4532	GL45	GL32	1
BD4538	GL45	GL38	1
BD3219	GL32	19/26	1
BD3224	GL32	24/29	1
BD3229	GL32	29/32	1




Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Dispensing cartridge

Model	Capacity	Materials	Pack Qty
BD1010	10ml	Piston PE, cylinder PP	3
BD1012	25ml	Piston PE, cylinder PP	3

Replacement parts

Model	Description	Pack Qty
BD1014	Discharge tube	1
BD1016	Filling tube, 250mm	2
BD1018	Filling tube, 500mm	2
BD1026	Seals	6



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.